

Sample Description:


Product Name : RP11304(SF-39)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP11304
 Lot No.: PJ00208506

	A	B
0.01min	20%	80%
25.0min	35%	65%
25.1min	100%	0%
30.0min		

STOP

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		19.763	500.960	4080.500	0.0253
2		20.905	929.527	12265.697	0.0762
3		21.775	5972.045	85135.117	0.5287
4		22.755	581251.563	15589398.000	96.8131
5		22.755	30196.592	353669.094	2.1964
6		26.093	2077.953	58027.398	0.3604
Total			620928.640	16102575.806	100.0000