


Sample Description:  
Summary:

Cat. No.: RP11324  
Lot No.: PX00207503

Product Name : RP11324(RR-8)  
Column : 4.6\*250mm, Hypersil-ODS2 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
0.01min A 18% B 82%  
25.0min 30% 70%  
25.1min 100% 0%  
30.0min STOP  
Flow rate : 1.0ml/min  
Wavelength : 220nm



### Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.000	584.430	28331.811	0.5095
2		15.027	1040.976	16024.151	0.2882
3		15.573	1003.869	15461.054	0.2780
4		16.148	418.712	4044.200	0.0727
5		17.003	140971.813	5356729.500	96.3309
6		17.003	9842.110	130154.969	2.3406
7		26.817	503.436	10015.396	0.1801
Total			154365.346	5560761.081	100.0000