


Summary:

Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient :
 A B
 0.1min 80 20
 25 min 55 45
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP11328

Lot No.: PL00104503



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		9.098	378.949	4420.994	0.2484
2		9.927	650.762	6232.232	0.3502
3		11.172	1827.335	23091.008	1.2974
4		11.528	145418.703	1746012.250	98.1040
TOTAL			148275.749	1779756.483	100.0000