


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 90 10
 25 min 0 100
 25.01min 0 100
 30 min 0 100

Cat. No.: RP11385
 Lot No.: PZ00220501

Flow rate : 1.0ml/min



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		7.562	385.217	2946.801	0.1164
2		14.597	455.075	2861.359	0.1131
3		15.367	880.318	6672.743	0.2637
4		17.640	345642.438	2518176.000	99.5068
TOTAL			347363.047	2530656.903	100.0000