


Sample Description:
Summary:

Product Name: RP11604 (AA-11)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A B
 11% 89%
25.0min 30% 70%
25.1min 100% 0%
30.0min
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP11604
Lot No.: PJ00521504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		16.047	368.370	2342.000	0.1602
2		16.220	1399.059	9336.600	0.6385
3		16.567	104823.836	1390996.500	95.1295
4		16.567	4274.109	19952.730	1.3646
5		17.550	3883.977	39585.504	2.7072
Total			114749.351	1462213.334	100.0000