


Sample Description:
Summary:

Product Name : RP11629(AC-AG-11)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 A B
 0.01min 7% 93%
 25.0min 18% 82%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat. No.: RP11629
 Lot No.: PJ00319505



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.538	4889.232	78799.648	0.4206
2		14.432	578653.250	18546256.000	98.9973
3		14.432	9925.424	20992.555	0.1121
4		16.227	5650.963	88064.898	0.4701
Total			599118.869	18734113.102	100.0000