


Sample Description:  
Summary:

Product Name: RP11738(TP-12)  
Column : 4.6×250mm, Sinochrom ODS-BP 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient

Cat. No.: RP11738  
Lot No.: PJ00525504

Flow rate : 1.0ml/min  
Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		16.295	460.150	2757.797	0.1702
2		16.697	137406.203	1546293.000	95.4264
3		18.203	4395.532	42197.105	2.6041
4		19.910	3086.000	29156.797	1.7994
Total			145347.885	1620404.700	100.0000