


Sample Description:

Summary:

Product Name : RP11739(LD12)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 A B
 0.01min 15% 85%
 25.0min 25% 75%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat No.: RP11739
 Lot No.: PJ00721504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		17.265	585.363	11858.499	0.0832
2		18.133	7500.079	154617.359	1.0847
3		19.200	306109.281	14087612.000	98.8321
Total			314194.723	14254087.858	100.0000