


Sample Description:

Product Name : RP11740(CK-12)
 Column : 4.6*250mm, Sinochrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient

Cat. No.: RP11740
 Lot No.: PJ00213505

	A	B
0.01min	15%	85%
25.0min	32%	68%
25.1min	100%	0%
30.0min	STOP	

Flow rate : 1.0ml/min
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.850	971.876	9905.205	0.3894
2		17.422	544.541	6201.604	0.2438
3		18.545	1914.342	22717.867	0.8930
4		18.832	151314.453	2485850.000	97.7187
5		19.987	441.459	4134.996	0.1625
6		22.362	1358.859	15074.995	0.5926
Total			156545.529	2543884.667	100.0000