


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 80 20
 25 min 55 45
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP12908
 Lot No.: PL00110503



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		11.618	1412.100	17078.398	0.3185
2		12.265	9911.313	113823.336	2.1224
3		12.902	398374.469	5149516.000	96.0200
4		13.582	5225.803	70375.539	1.3123
5		14.613	1017.581	12171.506	0.2270
TOTAL			415941.267	5362964.779	100.0000