


Sample Description:
Summary:

Product Name: RP12914(PEL-13-NH2)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A B
 14% 86%
25.0min 35% 65%
25.1min 100% 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP12914
Lot No.: PL00101502



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		16.265	1804.486	25939.689	1.4263
2		16.835	128411.266	1729189.250	95.0806
3		17.442	825.697	6691.900	0.3680
4		18.038	1189.122	6386.020	0.3511
5		18.168	3672.209	34314.875	1.8868
6		18.388	239.434	948.600	0.0522
7		18.567	1804.065	15185.704	0.8350
Total			137946.279	1818656.039	100.0000