


Sample Description:  
 Product Name : RP13314(IL-14-1-NH2)  
 Column : 4.6\*250mm, Kromasil C18 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient  
                   A                  B  
 0.01min      32%              68%  
 25.0min      45%              55%  
 25.1min      100%            0%  
 30.0min                        STOP  
 Flow rate : 1.0ml/min  
 Wavelength : 220nm

Cat. No.: RP13314  
 Lot No.: PL00301502



### Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.885	800.000	8857.899	0.2745
2		14.940	268.054	2629.604	0.0815
3		15.292	313.506	2755.100	0.0854
4		15.553	245.111	1850.101	0.0573
5		15.830	188.268	1498.098	0.0464
6		16.372	220.791	1976.147	0.0612
7		17.083	2614.462	14605.933	0.4526
8		17.295	151081.031	3193264.250	98.9412
Total			155731.224	3227437.131	100.0000