


Sample Description:

Summary:

Product Name: RP13326 (PYF10-NH2)
 Column : 4.6×250mm, Sinochrom ODS-BP 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 0.01min A B
 25.0min 25% 75%
 25.1min 40% 60%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

Cat. No.: RP13326
 Lot No.: PJ00629504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		12.747	1435.433	15205.100	0.6536
2		13.215	166563.078	2260429.500	97.1620
3		13.215	5300.086	50819.023	2.1844
Total			173298.598	2326453.623	100.0000