


Sample Description:

Summary:

Product Name : RP13686(MF-12-NH2)  
 Column : 4.6\*250mm, Kromasil C18 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient  
                   A                  B  
 0.01min      22%                  78%  
 25.0min      42%                  58%  
 25.1min      100%                 0%  
 30.0min                          STOP  
 Flow rate : 1.0ml/min  
 Wavelength : 220nm

Cat. No.: RP13686

Lot No.: PJ01225504



Results

| Peak No. | Peak ID | Ret Time | Height     | Area        | Conc.    |
|----------|---------|----------|------------|-------------|----------|
| 1        |         | 14.533   | 352.317    | 3361.768    | 0.0478   |
| 2        |         | 14.848   | 13304.853  | 228751.578  | 3.2543   |
| 3        |         | 15.437   | 1901.390   | 46148.637   | 0.6565   |
| 4        |         | 17.045   | 324308.281 | 6751020.500 | 96.0414  |
| Total    |         |          | 339866.841 | 7029282.483 | 100.0000 |