

PRODUCT INFORMATION

Version: 09/08/2008

dATP (100 mM) ≥99% HPLC**Cat. No.:** C01577**Name:** dATP (100 mM) ≥99% HPLC**Size:** 1 ml, 10 ml, 50 ml and 250 ml**Cas:** 1927-31-7**Synonyms:** 2'-Deoxyadenosine-5'-Triphosphate,
Sodium Salt**Molecular Formula:** $C_{10}H_{13}N_5Na_3O_{12}P_3$ **Molecular Weight:** 557.2**Storage:** Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.