

Taq Supermix

Cat. No.: E00033**Size:** 100 rxns**Version:** 08/03/2006

PRODUCT INFORMATION

Description

Taq Supermix is ready-to-use mixture of DNA Polymerase, salts, magnesium, and dNTPs for efficient amplification of templates in reactions. All you have to do is to add template and primers, and this will cut your reaction setup time by half. The *Taq* Supermix (Cat. No. E00033) comes as three vials, and each vial contains about 1.5 ml of solution, which is 100 PCR reactions (45 µl each).

Taq Supermix contents

22 mM Tris-HCl (pH 9.0), 11 mM KCl, 11 mM (NH₄)₂SO₄, 2.2 mM MgSO₄, 220 µM dNTP Mix, 0.11% Triton X-100, 50 U recombinant *Taq* DNA Polymerase/ml and stabilizer.

Storage

Store at -20 °C after receiving.

** For non-clinical research use only. **

Tel: 732-885-9188 Fax: 732-210-0262 E-Mail: info@genscript.com www.genscript.com