

PRODUCT INFORMATION**dNTP Mixture, 10 mM each**

Version: 10/11/2010

Cat. No. D0056**Name:** dNTP Mixture, 10 mM each**Size:** 0.5 ml

Description: 10 mM dNTP Mixture is suitable for use in PCR, sequencing, nick translation, fill-in, cDNA synthesis, TdT tailing reactions, and dilution of radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 10 mM dNTP Mixture

contains all four of the dNTPs, each at a final concentration of 10 mM.

Storage: Store at -20°C.

Quality control: By HPLC, purity >99.0%. Guaranteed to be free of nuclease (DNase, RNase) and phosphatases activity.