

PRODUCT INFORMATION**dNTP (25 mM each) ≥99% HPLC**

Version: 10/11/2010

Cat. No.: C01581**Name:** dNTP (25 mM each) ≥99% HPLC**Size:** 1 ml, 10 ml, 50 ml and 250 ml**Description:** 25 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-in, cDNA synthesis, TdT tailing reactions, and dilution of

radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 25 mM dNTP Mix contains all four of the dNTPs, each at a final concentration of 25 mM.

Storage: Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.