

PRODUCT INFORMATION**dNTP (10 mM each) \geq 99% HPLC**

Version: 09/08/2008

Cat. No.: C01582**Name:** dNTP (10 mM each) \geq 99% HPLC**Size:** 1 ml, 10 ml, 50 ml and 250 ml**Description:** 10 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-ins, cDNA synthesis, TdT tailing reactions, and the dilution of

radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 10 mM dNTP Mix contains all four dNTPs, each at a final concentration of 10 mM.

Storage: Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.