

CERTIFICATE OF ANALYSIS

Product	Urocortin III, human
Formula	$C_{185}H_{307}N_{53}O_{50}S_2$
M.W.	4137.93
Product Code	RP10100
Storage	Store at -20°C

Test Items	Specifications	Results
Apperance	White to off white powder	Conforms
Purity (HPLC)	> 95%	> 95%
Solubility Test	Insoluble in water. Dissolve peptide with 100-200 ul of DMSO then dilute up with desired buffer.	Complies

Caution:
For laboratory or further manufacturing use only.
Not intended for household use.

Laboratory Control