

PRODUCT INFORMATION

Rink Amide AM resin

Version: 10/11/2010

Cat. No.: C01573

Name: Rink Amide AM resin

Size: 25 g and 100 g

Description: The polymer matrix is copoly (styrene-1 % DVB), 200-400 mesh.

Synonyms:

4-(2',4'-Dimethoxyphenyl-Fmoc-aminomethyl)-phenoxy

acetamido-norleucylaminomethyl resin

Loading: 0.40 - 0.80 mmole/g resin, as determined photometrically from the amount of Fmoc chromophore liberated upon treatment of the resin with DBU/DMF

Storage: Store at room temperature without light and moisture

