

PRODUCT INFORMATION**Rink Amide AM resin**

Version: 06/30/2008

Cat. No.: C01573**Name:** Rink Amide AM resin**Size:** 25 g and 100 g**Description:** The polymer matrix is copoly (styrene-1 % DVB), 200-400 mesh.**Synonyms:**

4-(2',4'-Dimethoxyphenyl-Fmoc-aminomethyl)-phenoxy

acetamido-norleucylaminomethyl resin

Loading: 0.40 - 0.80 mmole/g resin, as determined photometrically from the amount of Fmoc chromophore liberated upon treatment of the resin with DBU/DMF**Storage:** Store at room temperature without light and moisture