

## PRODUCT INFORMATION

### Fmoc-L-Arg(pbf)-OH

Version: 07/10/2008

**Cat. No.:** C01546

**Name:** Fmoc-L-Arg(pbf)-OH

**Size:** 25 g, 100 g and 500 g

**Cas:** 154445-77-9

**Synonyms:**

N- $\alpha$ -Fmoc-NG-2,2,4,6,7-pentamethyldihydrobenzofuran-5-sulfonyl-L-arginine

**Molecular Formula:** C<sub>34</sub>H<sub>40</sub>N<sub>4</sub>O<sub>7</sub>S

**Molecular Weight:** 648.8

**Storage:** Store at room temperature without light and moisture

**Specification:**

Assay:  $\geq$  98.00% (HPLC)

Solubility: 1 mmole in 2 ml DMF clearly soluble

