

PRODUCT INFORMATION

Fmoc-L-Arg(pbf)-OH

Version: 07/10/2008

Cat. No.: C01546

Name: Fmoc-L-Arg(pbf)-OH

Size: 25 g, 100 g and 500 g

Cas: 154445-77-9

Synonyms:

N- α -Fmoc-NG-2,2,4,6,7-pentamethyldihydrobenzofura
n-5-sulfonyl-L-arginine

Molecular Formula: C₃₄H₄₀N₄O₇S

Molecular Weight: 648.8

Storage: Store at room temperature without light and moisture

Specification:

Assay: \geq 98.00% (HPLC)

Solubility: 1 mmole in 2 ml DMF clearly soluble

