

## PRODUCT INFORMATION

### Fmoc-L-Met-OH

Version: 10/11/2010

**Cat. No.:** C01556

**Name:** Fmoc-L-Met-OH

**Size:** 100 g and 500 g

**Cas:** 71989-28-1

**Synonyms:** N-a-Fmoc-L-methionine

**Molecular Formula:** C<sub>20</sub>H<sub>21</sub>NO<sub>4</sub>S

**Molecular Weight:** 371.5

**Storage:** Store at room temperature without light and moisture

**Specification:**

Assay: ≥ 98.00% (HPLC)

Solubility: 1 mmole in 2 ml DMF, clearly soluble

