

## PRODUCT INFORMATION

### Fmoc-L-Lys(Boc)-OH

Version: 10/11/2010

**Cat. No.:** C01555

**Name:** Fmoc-L-Lys(Boc)-OH

**Size:** 25 g, 100 g and 500 g

**Cas:** 71989-26-9

**Synonyms:** N-a-Fmoc-N-e-t.-Boc-L-lysine

**Molecular Formula:** C<sub>26</sub>H<sub>32</sub>N<sub>2</sub>O<sub>6</sub>

**Molecular Weight:** 468.5

**Storage:** Store at room temperature without light and moisture

**Specification:**

Assay: ≥ 98.00% (HPLC)

Solubility: 1 mmole in 2 ml DMF, clearly soluble

