

PRODUCT INFORMATION

BOP-Reagent

Version: 07/17/2008

Cat. No.: C01519

Name: BOP-Reagent

Size: 25 g, 100 g and 500 g

Cas: 56602-33-6

Synonyms:

Benzotriazole-1-yl-oxy-tris-(dimethylamino)-phosphonium

Molecular Formula: C₁₂H₂₂F₆N₆OP₂

Molecular Weight: 442.5

Appearance: White crystalline powder

Storage: Store at room temperature without light and moisture

Specification:

Assay: ≥ 99.00% min (HPLC)

Melting Point: 136-140°C

