


Positive

MS Spectrum

Intensity



Sample Information

User :xu
Date Acquired : 2008-8-1 10:32:43
Sample Name : [Nle4,D-Phe7]- α -Melanocyte Stimulating Hormone (MSH),amide
Cat No. : RP10658
Injection Volum : 1.000 μ L
Theoretical MW : 1647.00
Observed MW : 1647.00

Interface :ESI
Nebulizing Gas Flow :1.50 L/min
CDL Temp : 250 $^{\circ}$ C
CDL Volt : -25.0V
Block Temp : 200 $^{\circ}$ C

Prerod Bias :+4.5kv
Detector :1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH