

PRODUCT INFORMATION

Fmoc-L-Ala-OH

Version: 10/10/2008

Cat. No.: C01545

Name: Fmoc-L-Ala-OH

Size: 100 g and 500 g

Cas: 35661-39-3

Synonyms: N-a-Fmoc-L-alanine

Molecular Formula: C₁₈H₁₇NO₄

Molecular Weight: 311.30

Storage: Store at room temperature without light and moisture

Specification:

Assay: ≥ 98.00% (HPLC)

Solubility: 1 mmole in 2 ml DMF, clearly soluble

