


Sample Description:

Product Name : RP10224[(DF)W-8-NH2-3]
Column : 4.6*250mm, Kromasil C18 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min A 28% B 72%
25.0min 33% 67%
25.1min 100% 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

RP10224
PL00205509



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		17.467	632.984	9096.201	0.3419
2		18.022	281.500	3476.100	0.1307
3		19.138	1581.364	28807.094	1.0829
4		21.535	116863.734	2618792.250	98.4445
Total			119359.584	2660171.645	100.0000