


Sample Description:

Product Name : RP11102(DN-21)
 Column : 4.6*250mm, Kromasil C18 5um
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
 Solvent B : 0.1% trifluoroacetic in 100% water
 Gradient
 A B
 0.01min 16% 84%
 25.0min 20% 80%
 25.1min 100% 0%
 30.0min STOP
 Flow rate : 1.0ml/min
 Wavelength : 220nm

RP11102
 PX00209508



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		7.998	2596.500	18476.498	0.3026
2		8.918	10579.113	169064.953	2.7688
3		15.305	4995.610	61279.402	1.0036
4		16.042	218063.094	5807885.500	95.1171
5		17.832	3169.880	49330.504	0.8079
Total			239404.197	6106036.857	100.0000