


Product Name: RP10355 (YY-36-NH2)  
Column : 4.6×250mm, Sinochrom ODS-BP 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
    A                    B  
0.01min          30%          70%  
25.0min          40%          60%  
25.1min          100%         0%  
30.0min                    STOP  
Flow rate : 1.0ml/min  
Wavelength : 220nm

RP10355

PX00709508



### Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		10.412	9032.290	97687.586	0.7239
2		10.412	3733.199	15445.851	0.1145
3		11.482	31069.789	411567.156	3.0497
4		12.330	539464.188	12970543.000	96.1120
<b>Total</b>			<b>583299.466</b>	<b>13495243.593</b>	<b>100.0000</b>