


MS Spectrum



434.25 [M+H]⁺

Sample Information Item
Date and Time : 2005-5-30 19:31:14
User : xie
Sample : RP10861 RS-4
Inj. Volume : 1.000
MW : 433.42

RP10861
PY00105509

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0V
CDL Temp: 250°C
Block Temp: 200°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN