


M/S Spectrum

Positive

100



Sample Information Item

Date and Time : 2005-3-3 19:40:54

User : xie

Sample : RP11133 RP-7-1

Inj. Volume : 1.000

MW : 756.87

RP11133

PX00405509

Probe: ESI  
Nebulizer Gas Flow: 1.5L/min  
CDL: -20.0v  
CDL Temp: 250°C  
Block Temp: 200°C

Probe bias:  
Detector:  
T. Flow:  
B. conc:

+4.5kv  
1.5kv  
0.2ml/min  
50%H2O/50%ACN

m/z