


MS Spectrum



Sample Information Item : 2005-5-12 11:04:24  
Date and Time : xie  
User : RP10224((DF)W-8-NH2-3)  
Sample : PL00205509  
Inj. Volume : 1.000  
MW : 1132.38

Probe: ESI  
Nebulizer Gas Flow: 1.5L/min  
CDL: -20.0v  
CDL Temp: 250 C  
Block Temp: 200 C

Probe bias: +4.5kv  
Detector: 1.5kv  
T. Flow: 0.2ml/min  
B. conc: 50%H2O/50%ACN