


MS Spectrum



Sample Information Item

Date and Time : 2005-4-28 13:37:12

User : FJS

Sample : RP11102(DN-21)

Inj. Volume : 1.000

MW : 2436.64

RP11102

PX00209508

Probe: ESI

Nebulizer Gas Flow: 1.5L/min

CDL: -20.0V

CDL Temp: 250°C

Block Temp: 200°C

Probe bias: +4.5kV

Detector: 1.5kV

T. Flow: 0.2ml/min

B. conc: 50%H2O/50%ACN