


MS Spectrum



Sample Information Item

Date and Time : 2005-4-18 16:10:07
User : FJS
Sample : RP10289(RI-14)
Inj. Volume : 1.000
MW : 1633.93

RP10289
PL00109508

Probe: Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C

ESI

Probe bias: +4.5kv
Detector: 1.5kv
I. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN