


MS Spectrum



Sample Information Item
Date and Time : 2005-6-15 21:48:51
User : FJS
Sample : RP10467
Inj. Volume : 1.000
MW : 4243.75

RP10467
PY00409508

Probe bias: +4.5kv
Detector: I.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C