

Multifunctional Magnetic Separation Rack

Description:

GenScript Multifunctional Magnetic Separation Rack (M00140) is designed to facilitate efficient magnetic separation with 1.5-ml, 15-ml, or 50-ml centrifuge tubes. Using strong neodymium rare earth magnets, GenScript Multifunctional Magnetic Separation Rack can help the user achieve rapid separation of magnetic particles from the suspension. This rack can hold eight 1.5-ml tubes, one 15-ml tube, and one 50-ml tube.

Storage and Stability:

Magnetic separation racks should be stored in a cool, dry environment. They can be cleaned with a damp cloth and mild detergent. Do not autoclave.

Safety:

Please handle with care. Direct contact with the magnet material should be avoided, especially by persons with pacemakers or implants. Remove and dispose of any loose ferrous material.

Keep all magnetic media, watches, and sensitive electronic devices away from the magnetic separation racks. Information stored in credit cards, tapes and disks could be erased when placed in a strong magnetic field.

