


MS Spectrum



Sample Information Item
Date and Time : 2005-5-5 19:41:02
User : FJS
Sample : RP10355(YY-36-NH2)
Inj. Volume : 1.000
MW : 4309.85

RP10355
PX00709508

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
Block Temp: 200 C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN