


MS Spectrum



Sample Information Item : 2005-4-27 15:40:18  
Date and Time : xie  
User : RP10320 RA-8  
Sample : 1.000  
Inj. Volume : 958.14  
MW : RP10320 PX00609508

Probe: ESI  
Nebulizer Gas Flow: 1.5L/min  
CDL: -20.0v  
CDL Temp: 250°C  
Block Temp: 200°C

Probe bias: +4.5kv  
Detector: 1.5kv  
T. Flow: 0.2ml/min  
B. conc: 50%H2O/50%ACN