

PRODUCT INFORMATION

Version: 11/16/2008

BKC**Cat. No.:** C01683**Name:** BKC**Size:** 100 g**Cas:** 8001-54-5

Description: BKC is used as active pharmaceutical or cosmetic ingredient (API), Disinfection/Preservation. Benzalkonium application many include disinfecting instruments and preservative drugs in low concentration form.

Synonyms: Paralkan; Parasterol;

Alkyl-benzyl-dimethylammonium chloride; Ammonyx; Barquat MB-50; Barquat MB-80; Benirol; Benzalkonium Chloride; Bradophen; BTC; Cequaryl; Drapolene; Dropolex; Enuclene; Germitol; Gesminol; Osuan; Quaternary ammonium compounds

Formula: C₂₂H₄₀ClN

Molecular Weight: 354

Appearance: White Powder

Storage: Room temperature.