

## PRODUCT INFORMATION

Version: 11/16/2008

### BKC

**Cat. No.:** C01683

**Name:** BKC

**Size:** 100 g

**Cas:** 8001-54-5

**Description:** BKC is used as active pharmaceutical or cosmetic ingredient (API), Disinfection/Preservation. Benzalkonium application many include disinfecting instruments and preservative drugs in low concentration form.

**Synonyms:** Paralkan; Parasterol;

Alkyl-benzyl-dimethylammonium chloride; Ammonyx; Barquat MB-50; Barquat MB-80; Benirol; Benzalkonium Chloride; Bradophen; BTC; Cequartyl; Drapolene; Dropolex; Enuclene; Germitol; Gesminol; Osuan; Quaternary ammonium compounds

**Formula:** C<sub>22</sub>H<sub>40</sub>ClN

**Molecular Weight:** 354

**Appearance:** White Powder

**Storage:** Room temperature.