

## PRODUCT INFORMATION

Version: 12/04/2008

### BPTAB

**Cat. No.:** C01671

**Name:** BPTAB

**Size:** 100 g, 500 g, and 20000 g

**Cas:** 3779-42-8

**Description:** This peptide is a phase transfer catalyst, widely used as a lysing reagent in pharmaceutical and clinical diagnostic kits. In biotechnology and electronics, it is also used as an emulsifier, germicide,

and disinfectant, antistatic agent.

**Synonyms:** 3-Bromopropyltrimethylammonium bromide

**Molecular Formula:** C<sub>6</sub>H<sub>15</sub>Br<sub>2</sub>N

**Molecular Weight:** 261

**Appearance:** White crystalline powder

**Storage:** Store the product at room temperature.

