

**PRODUCT INFORMATION**

Version: 12/04/2008

**BPTAB****Cat. No.:** C01671**Name:** BPTAB**Size:** 100 g, 500 g, and 20000 g**Cas:** 3779-42-8

**Description:** This peptide is a phase transfer catalyst, widely used as a lysing reagent in pharmaceutical and clinical diagnostic kits. In biotechnology and electronics, it is also used as an emulsifier, germicide,

and disinfectant, antistatic agent.

**Synonyms:** 3-Bromopropyltrimethylammonium bromide**Molecular Formula:** C<sub>6</sub>H<sub>15</sub>Br<sub>2</sub>N**Molecular Weight:** 261**Appearance:** White crystalline powder**Storage:** Store the product at room temperature.