


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 02/23/2009 10:40:12
Sample Name : Beta-Amyloid A4 Protein Precursor(740-770),APP(C31)
Cat No. : RP20165
Lot No. : CB01903
Injection Volum : 1.000µL
Theoretical MW : 3717.10
Observed MW : 3717.00

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH