

CERTIFICATE OF ANALYSIS

Product information

Product name	Protein G
Cat. No.	Z02007
Lot No.	P01120601
Storage condition	Store at 4°C

QC results

Test Items	Specifications	Results
Purity	The purity is >90% by SDS-PAGE.	Pass
Stability	Place Protein G at 37°C for 7 days. The purity is >90% by SDS-PAGE.	Pass
Activity	Detect the activity using ELISA.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Richard* Date: 04/23/2009
Quality Assurance Department