

CERTIFICATE OF ANALYSIS

Product information

Product name	Protein L Resin
Cat. No.	L00239
Lot No.	R311904
Storage condition	Store at 4°C

QC results

Test Items	Specifications	Results
Binding capacity	≥ 15 mg/ml.	Pass
Stability	No decrease in binding capacity after incubated at 37°C for 1 week.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Richard*
Quality Assurance Department

Date: 04/03/2009