

CERTIFICATE OF ANALYSIS

Product information

Product name	High-Affinity Iodoacetyl Resin
Cat. No.	L00403
Lot No.	R311905
Storage condition	Store at 4°C

QC results

Test Items	Specifications	Results
Binding capacity	$\geq 3\text{mg/ml}$	Pass
Stability	Binding capacity keep above 3mg/ml after incubate at 37°C for 1 weeks.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Richard* Date: 05/08/2009
Quality Assurance Department