

CERTIFICATE OF ANALYSIS

Product information

| | |
|-------------------|-----------------|
| Product name | Protein G Resin |
| Cat. No. | L00209 |
| Lot. No. | R313904 |
| Storage condition | Store at 4°C |

QC results

| Test Items | Specifications | Results |
|------------------|---|---------|
| Binding capacity | $\geq 20\text{mg/ml}$ | Pass |
| Stability | Binding capacity keep above 20mg/ml after incubate at 37°C for 2 weeks. | Pass |

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Richard* Date: 04/13/2009
Quality Assurance Department